

Intubation Aid for Mice and Rats (small model)



Small model

This model is suitable for mice and small rats (\pm < 350gr)

How does it Work?

- Expose animal to anesthetics in induction chamber for initialanesthesia
- Take animal out of induction chamber and place "cord" around its frontteeth
- Slowly pull the animal with its snout/nose into the nosecone
- Secure code2utf('39',0)cordcode2utf('39',0)
- Because the intubation Aid is connected to the the animal is exposed toanesthetic gas constantly
- The "excessive /exhausted gas" is being disposed off through the connected Gas ExhaustUnit
- After intubation the animal can be connected to the UNO Micro

Ordernumber

Ordernumber	180000014
	Intubation Aid for mice and small rats $(\pm < 350gr)$







Use

- 1. Top part of the intubation stand with tread through "nose-cone".
- 2. Mouse placed with snout in nose-cone, Mousehan-ging by its front teeth on the tread supported by the back support.
- 3. Mouse hanging on intubation stand with anesthetic gas supply- and gas exhaust tubes connected to the intubationstand.
- 4. Checking the length of the intubation tube (length should be just before the bifurcation of thetrachea).
- 5. Light source placed at animal below the vocal cords of the animal. (Do not use a too strong light source because of to much light spreading).
- 6. Pull the animals tongue a little bit thus enabling you to see the vocal cords (Light comingthrough).
- 7. Place the intubation tube into the trachea in between the vocalcords.
- 8. Correct position of the intubation tube with regards to length of the tube and position of bifurcation in trachea.

The animal is now ready to be connected to the ventilator. Make sure that the tube maintains its position in the animal!