



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U102 AP-Microdriver (10μm)



This microdriver can be fine-tuned in the range of 25mm with 10μm steps. It will be even more accurate to work with the manipulator on the stereotax.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be