




CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U096 Device for Mice Spinal Card Fixation



This spinal adaptor can apply for mouse study. During the experiment, the mice can be lifted to avoid breathing effect. This adaptor can be rotated with different angles, depending on the height of the base, to facilitate observation under a microscope. Also, the base plate can be easily fixed to the microscope stage by T-shape screws.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be