



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U411 A/P Extension Assembly



A/P slide attachment can be easily installed on the operating arm of the V-shaped groove, and the sliding range is 60mm, with 0,1mm accuracy. When the electrode of probe is held perpendicular to the Z-axis, it can be manipulated in the horizontal direction, which provides fine control for targeting the coordinated of the brain atlas location.

(2-dimensional Manipulator not included)



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be