



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U399 Electrode Zero Point Calibrator



This electrode position calibrator is used to set the Zero Position for electrode in stereotaxic surgery. In the vertical direction, the reference point is 30mm higher than that of ear bars. It was designed for 68U000 stereotaxes.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be