

Cahteters, PC Series



Features

Pigtail connections

The PigTail connector features a titanium barbed connector pin attached to 2,5cm of Tygon tubing, a female luer lock, and an injection cap. The PigTail is inserted into the cathe- ter when an infusion is necessary.

Ordering Information

PC Series - PigTail Connectors	
PC-3	use with 3 Fr. catheter
PC-3,5	use with 3,5 Fr. catheter
PC-5	use with 5 Fr. catheter
PC-7	use with 7 Fr. catheter
PC-9	use with 9 Fr. catheter