



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U073 Head Rotation Platform for Rat & Mouse



The height of the platform can be adjusted according to the specific experimental requirements; the head holder can be adjusted in the vertical direction and 360° rotation horizontally. The platform below the head is suspended, especially for rats and mice that require 180° rotation in the head and body during surgery, such as experimental study on cerebral ischemia by thread metabolism, real-time monitoring of changes in the blood flow in the head (fixed fiber on the skull surface) to determine whether the model is successful or not.

Ordering number

68U071	Mouse rotation platform, with mouse rotation adaptor and 60° ear bar
68U074	Rat rotation platform, with rat rotation adaptor and 18° ear bar
68U073	Rat/Mouse rotation platform, with rat and mouse rotation adaptor and ear bar



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be