



Rodent Table for Rodent Surgery



Features

The non-porous, stain resistant, autoclavable, and durable Corian working surface of the rodent surgical table measures 20.3×29.2 cm, is 1.3cm thick and weighs 1.5kgs.

Grooves radiate from a central well for the drainage or placement of the flushing line during surgery.

A right angle Huber point flushing line enables flushing of the catheter and monitoring of catheter patency during surgery.

Four adjustable surgical rubber restraining straps and numerous hold-down slots along the perimeter of the table allow for the accommodation of various size rodents.