




**CARFIL QUALITY**  
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

## Taper Mask



This Taper Mask can easily be fastened on the platform and just place the animals nose into the mask hole during the experiment. The animal's position can be varied flexibly, such as lie down or lie flat without any influence on the research. Available in sizes to fit both mouse and rat.

Order No.	Product information
68U614	Rat Taper Mask
68U618	Mouse Taper Mask



**PAVAN SERVICE BVBA**

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

[www.carfil.be](http://www.carfil.be)