



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U078 Rat Anesthesia Adaptor (Tubing included)



The adaptor holds the rat head by mask and incisor bar, with table and good effect. The mask outlet can be connected to the gas exhaust unit to evacuate the non used and exhaled gas thus ensuring the safety of operators and environment. The connection between the head and the fixed plate ensures convenient vertical adjustment.

Note: Choose your required ear bars.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be